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CE 0197 **IP22**

1 INTENDED PURPOSE

The WBPM-1061 is a fully automatic, digital sphygmomanometer for monitoring blood pressure in the wrist. It enables quick and accurate measurement of systolic and diastolic blood pressure, as well as measurement of the pulse by means of the oscillometric method.

Please note that the results of the measurement are for your reference only, and they are not a substitute for diagnosis or treatment. If your results are abnormal, always seek the advice of a physician or medical professional.

The product is intended only for private household use by adults.

2 FUNCTIONALITY

- Automatic arm cuff inflation
- Adjustable cuff: 14 cm – 19.5 cm
- Memory storage for 2 users, each holding up to 60 measurements with date and time stamps
- Clear LCD with legible symbols
- Two AAA 1.5 V \equiv batteries (included) for power supply

3 TECHNICAL SPECIFICATION

| | |
|--|--|
| Model | WBPM-1061 |
| Lot No. | POJM210061 |
| Power supply | 3 V \equiv (2 x AAA 1,5 V \equiv battery (included)) |
| Measuring method | Oscillometric |
| Measuring range | Oscillometric Cuff pressure: 0-300 mmHg Pulse: 40-180 bpm Systolic pressure: 60-260 mmHg Diastolic pressure: 40-199 mmHg |
| Measuring accuracy (depending on temperature and humidity) | Blood pressure: \pm 3 mmHg Pulse rate tolerance: \geq 60 \pm 5% $<$ 60 \pm 3% |
| Arm cuff inflation | Automatic inflation with air (pump) |
| Ambient environment – operation | Temperature: from +10°C to +40°C Relative humidity: \leq 85%RH Atmospheric pressure: 80 kPa - 105 kPa (Measurement will be inaccurate under other conditions) |
| Ambient environment – transport and storage | Temperature: from -20°C to +50°C Relative humidity: \leq 85%RH |
| Device dimensions | 83mm x 74mm x 26mm |
| Cuff circumference | 14cm - 19.5cm |
| Application part type | BF |
| Operation mode | Continuous |
| Battery lifetime | Approx. 200 measurements |

4 SAFETY INSTRUCTIONS


1. Read the entire manual before use.
2. The blood pressure monitor is intended for adults and the patient should use the device himself.
3. The device is not intended for use on newborns, children and pregnant women. (No clinical trials have been conducted in newborns, children, or pregnant women.)
4. Consult a physician or other health care professional before older children use the blood pressure monitor.
5. The device cannot be used by people with severe cardiac arrhythmia.
6. Persons suffering from persistent arrhythmias, vasoconstriction, advanced atherosclerosis or diabetes mellitus and users of pacemakers should consult a physician before self-measurement of blood pressure, as in these cases variations in blood pressure may occur.
7. Swallowing batteries and / or battery fluid can be extremely dangerous. The device and batteries should be kept out of the reach of children and disabled

8. Strong magnetic fields, such as a cell phone, microwave oven, etc., should be avoided, as they may interfere with the operation of the device and adversely affect the measurement results.
9. It is recommended to keep the blood pressure monitor at least 30 cm away from other wireless devices such as a WLAN module, microwave oven etc. It cannot be used in the vicinity of an active HF SURGICAL EQUIPMENT and a ME SYSTEM RF shielded room. for magnetic resonance imaging, where EM INTERFERENCE intensity is high.
10. For information on potential electromagnetic or other interference between the monitor and other devices, see ELECTROMAGNETIC COMPATIBILITY INFORMATION, along with advice on how to avoid such interference.
11. Do not disassemble or repair the device. Always have the repair done by a professional or service technician.
12. Do not use a cuff other than the one provided by the manufacturer, otherwise it may cause a biocompatible hazard and cause a measurement error.
13. The use of this device for blood pressure measurement is not a substitute for a medical examination.
14. On the basis of the measurements made, you should not change the doses of medications prescribed by your doctor.
15. If you experience abnormal results, consult your doctor.
16. In the event of potential allergic reactions, stop using the device.
17. Do not clean, maintain, or repair the equipment while the pressure is being measured or turned on.
18. Do not use the device while driving any vehicle.
19. Keep this manual - it should be available to other users.
20. The electronic version of the manual is available at: <http://instrukcje.ver-shold.com>

5 EXPLANATION OF SYMBOLS

| | |
|----------------|--|
| | Symbol for warning |
| IP22 | The degree of protection of the electrical device against the penetration of external factors. The designation consists of the letters IP and two characters, the first of which indicates resistance to penetration by solids, and the second to water penetration. 2 protection against bodies larger than 12.5 mm (accidental finger contact), 2 protection against water dripping from various directions at angles up to 15 degrees |
| | Follow the operating instruction. |
| | Disposal of waste electrical and electronic equipment – see the DISPOSAL section of this manual. |
| | Instructions for the disposal of batteries – see the DISPOSAL section |
| | Designation of the packaging material – corrugated cardboard. |
| | Symbol indicating that the manufacturer has contributed financially to the construction and operation of packaging material recovery and recycling system. |
| SN | Serial number. |
| LOT | The symbol marks the lot number. |
| EC REP | The symbol marks the data pertaining to the authorised representative in Europe. |
| | The symbol marks the data pertaining to the Manufacturer. |
| | The symbol regarding the device year of manufacture. |
| | Means part of the BF class going into direct contact with the patient's body, meeting the requirements of the standard IEC 60601-1. The cuff is BF applied part. |
| CE 0197 | This product complies with the relevant EU directives. The number next to the "CE" mark is the number of the notified body which validated the manufacturer's quality assurance, manufacturing and final control systems for this product. |
| | Segregate the waste – Marking indicating the waste segregation necessity. |
| | Packaging: paper – Designation of the relevant waste bin. |
| | Packaging: plastic – Designation of the relevant waste bin. |

6 OVERVIEW

| Category | Systolic pressure (mmHg) | Diastolic pressure (mmHg) |
|-----------------------------------|--------------------------|---------------------------|
| Hypertension – stage 3 - severe | \geq 180 | \geq 110 |
| Hypertension – stage 2 - moderate | 160-179 | 100-109 |
| Hypertension – stage 1 - light | 140-159 | 90-99 |
| Above normal | 130-139 | 85-89 |
| Normal pressure | 120-129 | 80-84 |
| Normal pressure | $<$ 120 | $<$ 80 |

- Hypertension – stage 3 - severe
- Hypertension – stage 2 - moderate
- Hypertension – stage 1 - light
- Normal pressure

7 SET CONTENTS

| QUANTITY | Item |
|----------|--|
| 1 | Automatic blood pressure monitor with cuff |
| 1 | User's guide |
| 2 | AAA 1.5 V \equiv batteries |

Open the packaging and carefully take out the device. Make sure the set is complete, and that its components are undamaged. Check that the plastic parts and the tube are not damaged. If any part is damaged or missing, contact the dealer and do not use the product. Keep the package or recycle it in accordance with the applicable regulations.

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8 INSTALLATION AND CHANGING BATTERIES

CAUTION very careful! Read this manual.

- a. Open a battery compartment cover at the right down corner on the side of the monitor.

- b. Take out the used batteries if any.
- c. Load two new "AAA" size batteries paying attention to the polarity marked inside the battery compartment.

- d. Close the battery compartment cover.


When LCD shows battery symbol , replace all batteries with new ones. Rechargeable batteries are not suitable for this monitor.

Caution!

- The device is designed to be used with AAA 1.5 V \equiv batteries.
- The batteries are included.
- Observe the polarity (+/-) when inserting new batteries.
- Use the type of batteries recommended for this device.
- Do not put used batteries in a rubbish bin, but in special containers for used batteries.
- Never use spent batteries in combination with new ones.
- Never expose batteries to direct heat sources, such as excessive sunlight, heaters, or fire.
- Never throw batteries into fire – they may explode.
- When the device is not used, always remove the batteries. This will prevent damage to the device in the event electrolyte leaks from the batteries.
- Remove spent batteries from the device immediately.
- Spent batteries may leak, which may damage the device.
- If you touch acid from a battery, rinse your hands under running water. If any acid gets into your eyes, contact a doctor. The acid contained in the battery may cause irritation or burns.

9 SETTINGS OF TIME AND DATE

1. Once you install the battery or turn off the monitor, it will enter Clock Mode, and LCD will display time and date. See picture2.
2. While the monitor is in Clock Mode, press the "START" and "MEM" button simultaneously for two seconds, then release the buttons, the monitor will be in Clock and Date Adjustment Mode.
3. In Clock and Date Adjustment Mode, the time format will blink at first. See picture2-1.
4. While time format is blinking, press the button "MEM" to change the time format. The default time format is 24h.
5. Press the button "START" repeatedly, the year (range is 2018~2099), month, day, hour and minute will blink in turn. See picture2- 2&2-3 & 2-4 & 2-5 & 2-6.
6. While the number is blinking, press the button "MEM" to increase the number, keep on pressing the button "MEM", the number will increase faster.
7. During adjusting clock and date, the monitor will go back to Clock Mode automatically when no button will be pressed within 60 seconds.
8. You can turn off the monitor by pressing "START" button when the minute is blinking, then the time and date is confirmed.

Note:

- 1.1 The clock format could be set by user.
- 1.2 Table 2-1 instructs the conversion relations between 24-hour format and 12-hour format.

Table 2-1

| 24-hour format | 12-hour format | 24-hour format | 12-hour format |
|----------------|----------------|----------------|----------------|
| 0:00 | 12:00 AM | 12:00 | 12:00 PM |
| 1:00 | 1:00 AM | 13:00 | 1:00 PM |
| 2:00 | 2:00 AM | 14:00 | 2:00 PM |
| 3:00 | 3:00 AM | 15:00 | 3:00 PM |
| 4:00 | 4:00 AM | 16:00 | 4:00 PM |
| 5:00 | 5:00 AM | 17:00 | 5:00 PM |

| | | | |
|-------|----------|-------|----------|
| 6:00 | 6:00 AM | 18:00 | 6:00 PM |
| 7:00 | 7:00 AM | 19:00 | 7:00 PM |
| 8:00 | 8:00 AM | 20:00 | 8:00 PM |
| 9:00 | 9:00 AM | 21:00 | 9:00 PM |
| 10:00 | 10:00 AM | 22:00 | 10:00 PM |
| 11:00 | 11:00 AM | 23:00 | 11:00 PM |

10 BEFORE MEASUREMENT

- Read this manual before measuring blood pressure.
- Do not take a measurement immediately after a large meal. Wait at least one hour to obtain a more accurate reading.
- Prior to taking the measurement, at least for 30 minutes before do not smoke or consume alcohol, do not exercise intensively or undertake any intensive activity.
- If you are tired or exhausted, the reading can be non-representative for your normal condition.
- Rest for 5-10 minutes before taking the measurement.
- Sit down for the measurement and do not move or talk. Do not cross your legs.
- Install the cuff on the wrist approx. 2 cm above the palm.
- Sit comfortably during the measurement. The hand on which the pressure is being measured should rest at the table or desk with wrist at the height of the heart.
- The readings constantly change during the measurement. Minor changes of ambient or internal conditions (eating, drinking, talk, excitement, accelerated breathing, climatic conditions, etc.) result in blood pressure variation.
- Always measure pressure in the same hand (typically left). The pressure difference between your left and right hand can be significant.
- It is recommended to take daily measurements at the same time of day.
- Wait at least 5 minutes before taking the next measurement, otherwise the results can be inaccurate.
- Selecting the correct cuff size is very important for accuracy. The cuff included with the device is intended for adults.
- Taking the measurement too often may distort the blood circulation, resulting in injuries.
- Using the device for controlling blood pressure level will not result in prolonged impaired blood circulation.
- The device may be damaged or show incorrect results if it is exposed to unfavorable atmospheric conditions, such as humidity or extreme temperatures.
- Consult your doctor if you have any doubts about using the device, specifically in the following contexts:
 1. inflammation or other condition affecting the arm
 2. use of the cuff on the arm on the mastectomy side
 3. use of the product in combination with other health monitoring devices
 4. use of the cuff and its pressure on a limb, in which arteriovenous treatment is performed, or valves are applied
 5. checking, e.g. by observation, whether the use of the device does not result in prolonged blood circulation impairment in the user.

If irregular heart rhythm is detected during the measurement, the irregular pulse symbol (♥) will be displayed. Various factors can cause the irregularity. If the irregular pulse symbol is displayed, it is recommended that you consult a doctor for detailed examination.

There are two situations, in which the irregular pulse symbol (♥) can be displayed:

- 1) The heart rate variability coefficient > 25%.
- 2) The difference between the previous measurement ≤ 0.14 s, and the number of such heartbeats exceeds 53% of the total heartbeat number.

11 MEASURING BLOOD PRESSURE

Tips for taking measurements in seated position

1. Sit comfortably with your feet resting flat on the floor.
2. Do not cross your legs.
3. Lay your hand in front of you on an even surface (e.g. a table, a bench), with the palm inner side up, so that the cuff center is at the height of the heart.
4. The back should lean on the backrest.
5. Put on the cuff in accordance with the description: **Putting the cuff on**

Putting the cuff on

1. Fasten the cuff around a bare wrist, 1-2 cm above the wrist joint, on the inner side of the wrist.
2. While sitting, lay the hand in front of you on a desk or a table, with the cuff tightened around the wrist, and with the palm raised so that the cuff center is at the height of the heart. If the cuff is located correctly, you can see the readouts on the LCD.
3. The cuff cannot be too tight nor too loose.

UWAGA:

1. Always take measurements on the same wrist.
2. Do not move your arm, body, or the monitor during measurement.
3. Do not apply the cuff to an arm affected by inflammation, an acute disease, infection or skin injuries.
4. Stay calm for 5 minutes before the measurement.
5. Keep the cuff clean. If the cuff becomes dirty, clean it using a wet soft cloth and mild detergent.
6. Do not remove the cuff from the monitor.
7. Clean the cuff after each 200 uses.

Taking the measurement


CAUTION: To obtain a correct measurement with only minimum errors, the arm, palm and fingers should not move during the measurement. Do not move or talk during the measurement.

CAUTION: Please consult a health care professional for interpretation of pressure measurements.

1. After applying the cuff and your body is in a comfortable position, press the "START" button. All display characters are shown for self-test. You can check the LCD display according to the picture 6-1.
2. Then the current memory user choice will blink. See picture 6-2. Press "MEM" button to change over to other user. Confirm your selection by pressing "START" button. The current user can also be confirmed automatically after 5 seconds with no operation.
3. If the monitor has stored results, the LCD will momentarily display the most recent one. See picture 6-3. If no result has been stored, zero will appear on LCD. See picture 6-4.
4. Then the monitor starts to seek zero pressure. See picture 6-5.
5. Then the monitor inflates the cuff until sufficient pressure has built up for a measurement.
6. Then the monitor slowly releases air from the cuff and carries out the measurement.
7. Finally the blood pressure and pulse rate will be calculated and displayed on the LCD screen.
8. The blood pressure classification indicator and Irregular heartbeat symbol (if any) will blink on the screen. See picture 6-6. The result will be automatically stored in the monitor.
9. After measurement, the monitor will turn off automatically after 1 minute of no operation. Alternatively, you can press the "START" button to turn off the monitor manually.
10. During measurement, you can press the "START" button to turn off the monitor manually.

Picture 6-1

Picture 6-2

Picture 6-3

Picture 6-4

Picture 6-5

Picture 6-6

Reviewing historical measurements in the device "memory bank"

1. In Clock Mode, press "MEM" button, the monitor will display sign of current group. The amount of results in current user memory zone will be displayed. See picture 7-1. Press "START" button to switch group, press "MEM" to confirm current group. The current bank can also be confirmed automatically after 5 seconds with no operation.
2. Then LCD will display the average value of all results in the current bank. See picture 7-2. If no result stored in the current user memory zone, LCD will display "0" for blood pressure and pulse rate. See picture 7-3.
3. Press "MEM" button, the LCD will display the average value of the last three results in the current user memory zone. See picture 7-4. If no result stored in the current user memory zone, LCD will display "0" for blood pressure and pulse rate. See picture 7-5.
4. Press "MEM" button again, the memory number in the current user memory zone will be displayed. See picture 7-6. If the monitor has no result stored in the current user memory zone, the LCD will display "0". See picture 7-7.
5. Then the most recent result will be displayed with date and time stamp. See picture 7-8. Irregular heartbeat symbol (if any) and blood pressure classifica-

tion indicator will blink at the same time. Press "MEM" button repeatedly to review the results measured previously. If the monitor has no result stored in the current user memory zone, the LCD will display "0" for blood pressure and pulse rate. See picture 7-9.

6. When displaying the stored results, the monitor will turn off automatically after 1 minutes of no operation. You can also press the button "START" to turn off the monitor manually.

Picture 7-1

Picture 7-2

Picture 7-3

Picture 7-4

Picture 7-5

Picture 7-6

Picture 7-7


Picture 7-8

Picture 7-9

Deleting historical measurements from the device "memory bank"

When any result is displaying (exclude average value displaying), keeping on pressing button "MEM" for three seconds, all results in the current bank will be deleted. See picture 8. Press the button "START", the monitor will be turned off.

Picture 8

Technical alarm description

The monitor will show 'Hi' or 'Lo', as a technical alarm on LCD, with no delay, if the determined blood pressure (systolic or diastolic) is outside the rated range specified in part SPECIFICACIONES. In this case, you should consult a physician or check if your operation violated the instructions.

The technical alarm condition (outside the rated range) is preset in the factory and cannot be adjusted or inactivated. This alarm condition is assigned as low priority according to IEC 60601-1-8.

The technical alarm is non-latching and need no reset. The signal displayed on LCD will disappear automatically after about 8 seconds.

12 TROUBLESHOOTING

| PROBLEM | CAUSE | SOLUTION |
|--|---|---|
| The battery symbol is shown on the display | Low battery | Replace all the batteries |
| The display shows "Er0" | Arm or device movement during the measurement | Take another measurement, making sure that you do not move your hand or the device |
| The display shows "Er1" | Systolic pressure not detected | |
| The display shows "Er2" | Diastolic pressure not detected | Put the cuff on correctly, check the cuff connection with the main unit, and retry. |
| The display shows "Er3" | Cuff – wrong position or loose | |
| The display shows "Er4" | Cuff – wrong position or loose | |

| | | |
|---|--|--|
| The display shows "Er5" | Pressure > 300 mmHg | Wait 5 minutes and take another reading. If the display still shows any of the messages, contact the dealer and return the device for service. |
| The display shows "Er6" | Pressure > 15 mmHg for more than 160 seconds | |
| The display shows "Er7" | Internal memory error | Take the batteries out, wait five minutes and insert them again. |
| The display shows "Er8" | Device parameter error | |
| The display shows "ErA" | Pressure sensor parameter error | Put the cuff on correctly. Assume a correct body position. Repeat the measurement, stay calm, do not move, do not talk. It is recommended not to use the device. |
| No reaction after turning the device on or installing the batteries | Malfunction or strong electromagnetic interference | |
| The results on the display are incorrect | The cuff is put on incorrectly. The body position is incorrect. The patient moved his/her body, performed nervous movements, talked, etc. The user may suffer from arrhythmia. | |

13 MORE ABOUT BLOOD PRESSURE

Blood pressure is the force of blood against the artery walls as it circulates through the body. This pressure changes during each heartbeat cycle. The highest pressure is referred to as SYSTOLIC PRESSURE (SYS letters on the display). The lowest pressure is referred to as DIASTOLIC PRESSURE (DIA letters on the display). The doctor needs both these values to be able to evaluate the patient's blood pressure. The pressure values are affected by many factors, such as physical activity, a feeling of anxiety or the time of day. It changes continuously during the day. Blood pressure grows rapidly in the early morning, drops before noon, rises again in the afternoon, and drops to a minimum for the night. It can also change over a short period. Accordingly, the results of a series of measurement can be different.

A measurement taken at a clinic provides only a momentary value (snapshot). On the other hand, regular measurements taken at home provide a realistic picture of the blood pressure in the daily living environment. In addition, pressure measured at home is often different than that measured at a clinic because the individual concerned is more relaxed and comfortable. Regular measurements taken at home give the doctor valuable information on the patient's blood pressure in normal conditions.

The blood pressure values are typically given in millimetres of mercury mmHg. It is a unit of pressure exerted at the base of a 1-millimetre high column of mercury at 0°C under normal earth gravity. The device does not support measurement in kilopascals (kPa), where 1 mmHg = 133.3224 Pa.

The following matrix developed by the World Health Organization is typically used for estimating the health condition based on pressure, without addressing the person's age.

| Category | Sysstolic pressure (mmHg) | Diastolic pressure (mmHg) |
|-----------------------------------|---------------------------|---------------------------|
| Hypertension – stage 3 - severe | ≥ 180 | ≥ 110 |
| Hypertension – stage 2 - moderate | 160-179 | 100-109 |
| Hypertension – stage 1 - light | 140-159 | 90-99 |
| Above normal | 130-139 | 85-89 |
| Normal pressure | 120-129 | 80-84 |
| Normal pressure | < 120 | < 80 |

14 CLEANING AND MAINTENANCE

Correct and regular cleaning will improve the safety of the product and extend its lifetime. Switch off the device prior to cleaning and maintenance.

- Clean the device with a dry soft cloth.
- Do not clean the device with any chemicals, alkali, abrasive or disinfecting agents, as they may damage the surface of the device.
- Keep the cuff clean. If the cuff is dirty. If the cuff gets dirty, clean it using a wet soft cloth and mild detergent. Do not remove the cuff from the monitor. Do not dry the cuff in a dryer and do not iron it. It is recommended to clean the cuff every 200 measurements.
- The cuff Velcro maintains its properties up to 1000 fastening/opening cycles.
- The device keeps its safety and high quality of operation up to 10,000 measurements or for 3 years, if it is used from time to time.
- Protect the device from direct sunlight, humidity, and dust.
- It is recommended to check the device's performance every 2 years, and after each repair. Contact the service to have the device performance checked.
- If required (e.g., in a hospital or clinic environment), disinfect the cuff twice a week. Wipe the inner side (the one contacting skin) with a soft cloth moistened with a few drops of ethyl alcohol (75-90%), and leave the cuff to dry.

Caution!

Remove the batteries before prolonged storage of the device. Otherwise, the device may be damaged from possible electrolyte leak.

15 REPAIR

The device does not contain any user-serviceable parts. Do not attempt to repair the device on your own. Have the service or a professional perform the repair. Circuit diagrams, part lists, descriptions, calibration instructions, or other information determined by the customer can be submitted to the professional or service company.

Warning: Any modification of the product is prohibited.

16 HANDLING AND STORAGE

- It is recommended that the device be stored in its original packaging.
- Always keep the device in a dry, well-ventilated location inaccessible to children.
- Protect the product against extreme temperatures, humidity, dust, and direct sunlight.
- Electrolyte leak from the batteries can damage the device. Remove the batteries before prolonged storage of the device.
- Remove the batteries before transporting the product.
- Protect the device against vibration and shock during transport.
- If the device was stored at a low temperature (0°C or less), wait until it accommodates to the room temperature before using it. **The monitor requires 6 hours to warm from the minimum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 20 °C.**

17 DISPOSAL

All the packaging materials are 100% recyclable, and are labelled as such. The packaging should be disposed of in accordance with local regulations. Keep the packaging materials out of children's reach as these materials can pose a hazard.

Proper disposal of the device:

- According to the WEEE Directive 2012/19/EU, the crossed-out wheelee bin symbol (shown on the side) is used to label all electric and electronic devices requiring segregation.
- Do not dispose of the spent product with domestic waste: dispose of it at an electric and electronic devices collection and recycling point. The crossed-out wheelee bin symbol placed on the product, instruction manual or package communicates this requirement.
- The plastics contained in the device can be recycled in accordance with the specific marking. Thanks to recycling, re-use of materials or other forms of utilizing used equipment you make a significant contribution to the protection of our natural environment.
- Information on electric and electronic devices collection locations is available from local government agencies or from the dealer.

- Spent or fully discharged single-use and rechargeable batteries must be discarded in dedicated labelled bins, turned over to hazardous waste collectors or returned to electric equipment dealers.

18 CE DECLARATION

This product has been designed, manufactured and marketed in accordance with the requirements of the medical device directive (MDD) 93/42/EEC, and for that reason the product has been marked with the CE mark, and a declaration of conformity has been issued for it, which is made available to market regulators. The product complies with the following European standards applicable to non-invasive blood pressure measuring devices:
 EN 1060-1: Non-invasive sphygmomanometers. Part 1: General requirements.
 EN 1060-3: Non-invasive sphygmomanometers. Part 3: Additional requirements for electro-mechanical blood measuring systems.
 EN 60601-1: Medical electric devices. Part 1: General requirements for basic safety and essential performance.
 EN 60601-1-2 - Medical electric devices. Part 1-2: General requirements for basic safety and essential performance. Collateral standard: Electromagnetic disturbances – Requirements and tests”
 EN 60601-1-11 - Medical electric devices. Part 1-11: General requirements for basic safety and essential performance. Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.
 EN 60601-1-6 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Supplementary standard: Serviceability
 EN 60601-1-8 Medical electrical equipment - Part 1-8: General requirements for basic safety and essential performance - Supplementary standard: General requirements, tests and guidelines for alarm systems in medical electrical equipment and medical electrical systems
 IEC 80601-2-30 Medical electrical equipment –Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers)
 ISO 81060-2 Non-Invasive Sphygmomanometers - Part 2: Clinical Validation Of Automated Measurement Type.

19 ELECTROMAGNETIC COMPATIBILITY STATEMENT

Table 1
Guidance and manufacturer's declaration – Electromagnetic compatibility: Immunity test

| Phenomenon | Standard | Electromagnetic environment |
|----------------------------------|---------------------------|------------------------------|
| Electromagnetic emission tests | CISPR 11 Group 1, Class B | Home health care environment |
| Harmonic distortion | IEC 61000-3-2 | Home health care environment |
| Voltage fluctuations and flicker | Class A | Home health care environment |

Table 2
Guidance and manufacturer's declaration – Electromagnetic compatibility: Immunity test

| Phenomenon | Standard | Damage resistance test |
|---|---------------|--|
| | | Home health care environment |
| Electrostatic discharge | IEC 61000-4-2 | ± 8 kV on contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV via air |
| Emitted high-frequency electromagnetic field | IEC 61000-4-3 | 10V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz |
| Radiated RF electromagnetic fields in close proximity | IEC 61000-4-3 | See table No. 3 |
| Rated power-frequency magnetic fields | IEC 61000-4-8 | 30A/m 50Hz or 60Hz |

Table 3
Guidance and manufacturer's declaration – Electromagnetic compatibility: Radiated RF electromagnetic fields in close proximity

| Frequency (MHz) | Frequency span (MHz) | Damage resistance level |
|-----------------|----------------------|--|
| | | Healthcare facility professional environment |
| 385 | 380-390 | Pulse modulation 18 Hz, 27 V/m |
| 450 | 430-470 | FM, deviation ±5 kHz, 1 kHz, 28 V/m |
| 710 | 704-787 | Pulse modulation 217Hz, 9V/m |
| 745 | | |
| 780 | | |
| 810 | 800-960 | Pulse modulation 18Hz, 28V/m |
| 870 | | |
| 930 | | |
| 1720 | 1700-1990 | Pulse modulation 217Hz, 28V/m |
| 1845 | | |
| 1970 | | |
| 2450 | 2400-2570 | Pulse modulation 217Hz, 28V/m |
| 5240 | 5100-5800 | Pulse modulation 217Hz, 9V/m |
| 5500 | | |
| 5785 | | |

20 WARRANTY

In order to submit a complaint about the product, submit the device to the Customer Service Point at any Biedronka store. Should you have any questions or issues related to product operation, send them to the following e-mail address:

- infolinia@vershold.com
- or contact us by phone: +48 667 090 903

Information on the repair status may be obtained by contacting the relevant service desk.

Warranty and post-warranty service is provided by QUADRA-NET Sp. z o. o.

- tel. 61/853 44 44
- tel. kom. 664 44 88 00
- email: infolinia@quadra-net.com
- www: www.quadra-net.pl

Your opinion is important to us. Evaluate our product at: www.vershold.com/opinie

General Terms and Conditions of the Warranty

- The Warrantor for this product warrants that within 36 months from the

product purchase date, confirmed with a duly completed warranty card, the product shall be, at the discretion of the Warrantor, subject to:

- repairing the product free of charge, subject to the condition that any dismantled, damaged or faulty parts, in order to ensure the Customer's safety, shall become the Warrantor's property,
 - replacing the product with a new one at the point of purchase, subject to the condition that the damaged or faulty product shall become the Warrantor's property,
 - returning the price of the product to the User, subject to the condition that the damaged or faulty product shall become the Warrantor's property.
- A defective/damaged product is a product failing to provide the features described in the User's Manual due to intrinsic properties of the device.
 - The Beneficiary under the warranty shall deliver the appliance to the Customer Service Point at a given store in order to submit a complaint about the equipment. When asserting your rights under warranty, you need to present the proof of purchase (purchase receipt, invoice) together with a description of the defect in as much detail as possible, in particular of the external manifestations of the defect.
 - A warranty card without the vendor's stamp and signature or sale date is invalid.
 - The warranty excludes products with physical damage other than that caused by the manufacturer or distributor, including specifically any damage resulting from incorrect usage or force majeure.
 - The warranty excludes products without the tamper-proof seal, unless an authorized Service employee has removed the seal.
 - The repair center shall make every effort to complete the repair within 21 days. The said time limit may be extended to 1 month in the event that spare parts currently unavailable at the Service must be purchased.
 - It is recommended that the complete product be submitted together with the complaint so as to facilitate service operations.
 - This warranty for the sold consumer product does not exclude, restrict or suspend any rights of the Buyer related to laws and regulations concerning warranty for defects in sold items pursuant to the provisions of the Civil Code of 23 April 1964 (Journal of Laws Dz.U.2014.121).
 - The warranty shall be applicable within the territory of the Republic of Poland.

Warranty card

- Stamp
- Sale date
- Signature
- Appliance name
- Appliance model
- Serial No. or lot No.

Manufacturer (Warrantor)
 Andon Health Co., Ltd.
 No. 3 Jinping Street,
 YaAn Road, Nankai District,
 Tianjin, 300190, PRC

AUTHORIZED REPRESENTATIVE IN EUROPE:

iHealthLabs Europe SAS
 36 rue de Ponthieu,
 75008, Paris, France

Importer:

Vershold Poland Sp. z o.o.
 ul. Żwirki i Wigury 16A,
 02-092 Warsaw, Poland
 Made in China