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
AUTOMATIC BRACHIAL BLOOD PRESSURE MONITOR WITH VOICE INDICATION

Table of Contents

1. INTRODUCTION.....	28
2. INTENDED USE.....	28
3. TECHNICAL DATA.....	29
4. PRINCIPLES OF SAFE USE.....	29
5. EXPLANATIONS OF THE SYMBOLS.....	33
6. THE CONTENTS OF THE SET.....	35
7. CONSTRUCTION.....	36
8. USE.....	37
9. CLEANING AND MAINTENANCE.....	45
10. TROUBLESHOOTING.....	46
11. STORAGE AND TRANSPORT.....	47
12. WARRANTY AND SERVICE.....	48
13. DECLARATION OF CONFORMITY WITH EUROPEAN STANDARDS.....	48
14. MANUFACTURER AND IMPORTER:.....	48
15. ELECTROMAGNETIC COMPATIBILITY INSTRUCTIONS.....	48

1. INTRODUCTION

Thank you for choosing our product. We are confident that its use will be a great pleasure for you. Please read the user manual carefully before using the product, in particular the instructions on safe use. Please make it available to everyone who is going to use the device. Please keep this manual in order to be able to refer to it while using the product.

The colours of the product may slightly differ from the product shown in the illustrations in the manual and on the packaging.

2. INTENDED USE

This device is intended for private use and may not be used for commercial purposes. The device is designed to measure blood pressure (systolic and diastolic) as well as heart rate. Memory capacity allows to save up to 192 test results for each of the two users. Voice messages allow for easier operation even for people with sight problems. The blood pressure measured by this product is equivalent to that measured by auscultation of an artery. The device can be used even when the pulse is irregular, i.e. it differs by 25% from the average detected during the measurement. If the instrument detects an irregular pulse and a blinking symbol appears on the display, it is recommended to see a doctor. This device is designed to monitor blood pressure and pulse rate of adults in a domestic environment.

What is blood pressure?

Arterial pressure is the pressure exerted by the flowing blood on the arterial walls. Systolic blood pressure occurs when the heart contracts and the blood is pumped into the veins. Diastolic pressure when the blood returns to the heart. Blood pressure is measured in manometric unit of pressure - millimetres of mercury (mmHg).

What is hypertension and how to control it?

Hypertension is an extremely high blood pressure. Blood pressure should be checked regularly to detect hypertension as early as possible, as disregarding this can lead to serious illnesses such as heart attacks and strokes. To prevent or reduce the hypertension:

- do not smoke
- avoid stress
- maintain correct weight
- exercise regularly
- reduce consumption of salt and fats

Why is it worth to measure blood pressure at home?

Pressure measurement in a hospital or doctor's office may cause anxiety (the so-called white coat syndrome) and affect the increase in pressure by up to 25-30 mmHg. Measurement at home reduces all external factors.

3. TECHNICAL DATA

Product	Automatic Brachial Blood Pressure Monitor With Voice Indication
Model	KF-65A
Work method	Automatic pump, automatic measurement
Measuring procedure	Oscillometric method
Display	LCD Digital display
Cuff pressure range	0 – 280 mmHg
Pulse range	40 – 199 beats/min
Cuff pressure range	0 – 280 mmHg
Over-pressure protection	280 mmHg
Measuring accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Pressure measurement	Pressure transducer
Content	Monitor unit, arm cuff, instruction manual
Power source	6 V \approx , 4 AA x 1.5 V = alkaline battery Or optional AC adapter with micro-USB cable - output 5V = 500 mA - not included
Monitor Unit Dimensions	140 mm * 110 mm * 80 mm
Weight	393 g (with arm cuff, excluding batteries)
Cuff circumference	23 – 33 cm
Memory	2*192
Operating temperature / humidity / air pressure	+5 - +40°C / 15% - 85%RH / 86 - 106 kPa
Storage temperature / humidity / air pressure	-20 - 55°C / 15% - 93%RH / 86 - 106 kPa
Applied part	 = Type BF
Protection against electric shock	Internally powered ME equipment (When using only the batteries) Class II ME equipment (AC adapter powered, adapter not included)
Expected service life	5 years
Software version	001

4. PRINCIPLES OF SAFE USE

Read the entire user manual carefully before use. Keep it for future reference.

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State, where the patient or user resides.

 **CAUTION!**

- Do not charge non-rechargeable batteries (other than accumulators).
- Do not short-circuit the battery power terminals.
- Use batteries of the same type as recommended for this device.
- When installing new batteries, remember to check the +/- polarity of the batteries.
- Remove discharged batteries immediately from the unit, as they may leak and cause damage.
- Always remove the batteries when not in use - this will prevent the product from being damaged if the electrolyte leaks from the battery.
- Always replace the entire battery pack with a functional one.
- Keep batteries away from children and animals.
- Do not dispose of your used batteries in a waste bin, but in special containers for used batteries and accumulators.
- Avoid exposing the battery to very high/low temperatures (below 0°C/32°F and over 40°C/104°F). Extreme temperatures can affect their service life.
- Avoid contact with liquids and metal objects, as this may result in total or partial damage to the battery.
- Place worn or damaged batteries in a special container.
- Never expose batteries to direct sunlight, radiators, fire or excessive heat.
- Never dispose of batteries into a fire - they may explode.
- Never immerse batteries in water. Do not touch the batteries with wet hands.
- If hands come into contact with acid from the battery, rinse hands under running water. If electrolyte gets into your eyes, contact your doctor. Acid contained in batteries may cause skin irritation or chemical burns.
- Swallowing batteries can be fatal! Keep batteries away from children and pets. Immediately seek medical attention if the battery is swallowed.

 **Warning!!** *In order to ensure adequate safety measures and operating conditions, the battery should be mounted only by adults.*

1. Always abide by any kinds of prohibitions, regulations and warnings conveyed by the staff of the entity where you wish to use the device.
2. Do not use the device near flammable fabrics.
3. This appliance may be used only by adults.
4. Children must not play with this device.
5. Unattended children must not clean or maintain the unit. Cleaning and maintenance of the appliance must be carried out as described in section: 9. CLEANING AND CONSERVATION.
6. Avoid exposing the unit to very low or very high temperatures (below 5°C / 41°F or above 40°C / 104°F during work and below -20°C / -7,6°F or above 55°C / 131°F during storage and transportation). Extreme temperatures can affect service life of the device.
7. The device should not be used in dusty conditions. Make sure that animals do not have access to the product. Damaged product can have cracks or broken parts with dangerous sharp edges. Do not expose the device to sunlight and do not place it in close proximity to heat-generating devices. Plastic parts may deform, which may affect the operation of the device or cause it to break down completely.
8. Do not immerse the device in water or other liquids, as wetting will harm the electronics and may damage the equipment.
9. Do not use the unit if it has fallen to the ground or is visibly damaged.
10. The manufacturer is not responsible for any problems caused by unauthorized modifications to the equipment.
11. Only qualified personnel or an authorized service centre may install or repair this product. Repairing by an unauthorized or unqualified service may cause device damage and warranty loss. The manufacturer is not responsible for problems caused by unauthorized modifications to the equipment made by the user.
12. The device is used to measure and control of blood pressure. It can be used solely on the upper-arm. It cannot be used on other parts of the body or for purposes other than blood pressure measurement.
13. Taking a measurement is not equivalent to a professional diagnosis. The diagnosis can only be made by a doctor or a trained and authorized person. Always consult your doctor if you are treating abnormal blood pressure or heart rate.
14. If you are taking medication, it is important to consult your doctor to determine the right time to take your blood pressure measurements.

15. You are not allowed to change the dosage of medicines by yourself, without consulting your doctor, only on the basis of the measurement.
16. If you suffer from abnormal heart rhythm (premature ventricular contractions, atrial contractions, atrial fibrillation), the device may give erroneous results. Be sure to consult your doctor about the results.
17. Make sure that the device is not exposed to magnetic field interference, as this may affect the measurement results.
18. It is not advisable for pregnant women or women suspecting of being pregnant to use the device. The results may be distorted and the influence of such devices on the foetus is not known.
19. The device should not be used for continuous measurements during procedures, as it may cause changes in pressure and swelling in the limb on which the cuff is placed. This can also lead to cardiac arrest.
20. The material used for the cuff is intended for direct contact with the skin and should not cause any sensitization or irritation. However, if any skin reactions occur, immediately stop using the cuff and device and consult a physician as soon as possible.
21. Use accessories and equipment from the set or received from the manufacturer (service/replacement), because otherwise erroneous measurements may occur, the user's health may be endangered or the device may be damaged.
22. Continuous cuff compression due to kinking of the connection tubing can lead to blood flow interference and subsequent patient injury.
23. Placing the cuff and compressing it on any limb where there is intravascular access or therapy or an arterio-venous (A-V) shunt may lead to temporary interference with blood flow and may result in patient injury.
24. Avoid pumping the cuff on the mastectomy side of the arm.
25. Inflating the cuff may cause temporary loss of function of concurrently used ME equipment placed on the same limb as the cuff.
26. The automatic blood pressure monitor needs to be checked to ensure that its operation does not cause prolonged impairment to the patient's blood circulation.
27. Small parts – choking hazard.

5. EXPLANATIONS OF THE SYMBOLS

	Warning!! To reduce the risk of injury, the user must read the user manual.
LOT	Lot number
SN 612736XXXXX	Device serial number.
	BF protection type. Protection against electric shock.
	Caution.
	Important information on the use of the device.
	Warning.
	CE0123 - 0123 denotes the number of the notified body examining the product.
	Keep dry.
	Keep away from sunlight
	This way up
	Handle with care
	Manufacturer
EC REP	Sales representative within the EU.

	<p>The packaging of the device can be recycled in accordance with local regulations.</p>
	<p>The packaging components must be disposed of in the appropriate containers</p>
	<p>Marking of the container to which the foil should be thrown - plastic / metal.</p>
	<p>Designation of the container to which the packaging should be thrown - paper.</p>
	<p>Symbol for cardboard (packaging material).</p>
	<p>In accordance with directive 2013/56/UE on the disposal of batteries and accumulators, this product is marked with the symbol of a crossed-out waste bin. The symbol indicates that the batteries used in this product should not be disposed with regular household waste, but treated according to the legal directives and local regulations. It is not allowed to dispose of batteries and accumulators with unsorted municipal waste. Users of batteries and accumulators must use the available collection points of these items, which allows them to return, recycle and disposal. Within the EU, the collection and recycling of batteries and accumulators is a subject to separate procedures. To learn more about existing procedures of recycling of batteries and accumulators please contact your local office or an institution for waste disposal or landfill.</p>
	<p>The device is marked with a crossed-out garbage bin, in accordance with the European Directive 2012/19/UE on used electrical and electronic Equipment (Waste Electrical and Electronic equipment - WEEE). Products marked with this symbol should not be disposed of or dumped with other household waste after a period of use. The user is obliged to get rid of used electrical and electronic equipment by delivering it to designated recycling point in which dangerous waste are being recycled.</p> <p>Collecting this type of waste in designated locations, and the actual process of their recovery contribute to the protection of natural resources. Proper disposal of used electrical and electronic equipment is beneficial to human health and</p>

	environment. To obtain information on where and how to dispose of used electrical and electronic equipment in an environmentally friendly manner user should contact their local government office, the waste collection point, or point of sale, where the equipment was purchased.
--	--

6. THE CONTENTS OF THE SET

After unpacking, check that all replaced parts are in the package and that the device is working properly. Included in the set:

- Automatic brachial blood pressure monitor with voice indication
- Cuff with an airway tube.
- Storage case for the device
- 4 x AA 1.5 V battery =
- User manual


 *If any of the above items is missing or damaged, please contact your dealer.*

 *We advise to keep the container in case of any complaints. Keep the container out of reach of children.*

You can obtain a free printed version of the manual with the large font size or download it from the web address below:

<https://mpotech.eu/qr/hoffen>


or by scanning the QR code.


7. CONSTRUCTION

APPEARANCE OF THE DEVICE AND ITS COMPONENTS

Blood pressure monitor base:

Cuff with air tube ending in a plug:


 All operations should be carried out gently so as not to damage the components.

 4 x 1.5 V AA battery included in the set =.

INSTALLATION AND REPLACEMENT OF BATTERIES

 Incorrectly inserting the battery will only result in a malfunctioning device.

Unpack the device. Check whether any of the components are damaged. Place the device on a soft and delicate surface (blanket, towel, etc.) so that you have free access to the battery compartment.

Open the battery compartment cover by sliding and removing the cover as illustrated. When replacing the batteries, remove the entire set of 4 batteries and dispose of them according to the regulations.

Insert the whole new set of 4 AA 1.5V batteries = according to the polarity marked with the drawings inside the chamber.

Close the battery compartment cover until the latch closes, preventing the batteries from falling out.

TAKING A MEASUREMENT

 Immediately before the pressure measurement, we recommend that you do not eat, drink, smoke or exercise. Sit down and relax for about 10 minutes before the test.

 Remove clothing tightly adjacent to the arm on which the cuff is to be worn. In order to maintain the same measurement conditions, to compare values, the cuff must be worn over the same arm. We recommend that you take measurements on your left shoulder at the same time of day, if possible.

 Prior to measurement, the cuff on the arm in which the muscles are relaxed should be carefully adjusted. Incorrectly adjusted and fastened cuffs may cause an incorrect measurement.


 If the measurement is to be repeated, wait 2-3 minutes between them. It is best to remove the cuff and lift the arm so that the amount of blood and its flow in the arm can normalize. The measurement results may be incorrect otherwise.

 If it is not possible to adjust the cuff on the left upper arm, it can be placed on the right upper arm as well. However, all measurements shall be made on the same arm.


 Measurement results may be influenced by the measurement site, patient position, exercise, or physiological condition of the patient.

 Too frequent measurements can cause injury to the patient due to blood flow interference. Placing the cuff over a wound may cause further injury.


PREPARATION AND APPLYING THE ARM CUFF

1. Insert and push the airway plug from the cuff gently into the device **shown in the illustration above**

2. If the cuff is already prepared, ignore this point. Insert the end of the cuff (including the rubber clip) into the metal band so that it forms a loop. The Velcro should be on the outside.
3. Place the cuff on the left shoulder so that the tube extends towards the forearm - see above.
4. Attach the cuff to the shoulder as shown in the illustration above.
5. Make sure the lower edge of the cuff is about 2-3 cm above the elbow and that the tube is on the inside of the arm.
6. Wrap the cuff firmly in place around your arm using the cloth fastener as shown in the previous illustration. Air tube coming out of the cuff should be located above the brachial artery (~3cm long strap), on the inside of the arm.
7. There should be no space between the cuff and arm, as this may affect the readings. Rolled up sleeves must not press down on the shoulder and must therefore be rolled up or removed. The adjustable cuff size allows it to be applied to an arm circumference of 23 to 33 cm.
8. Secure the cuff with Velcro so that it is comfortable and not too tight. Remember to sit comfortably with your legs uncrossed and feet flat on the floor, back and arm supported. Lay your arm on the table (palm up) with the centre of the cuff on the same height as the right atrium of your heart as shown in the figure below.

9. Make sure that the air tube is not bent or pressed.
10. Stay in this position for five minutes before starting the measurement.

TAKING A MEASUREMENT

Once the cuff is correctly positioned on the arm, the pressure measurement can be started. After the correct measurement, the result will be stored automatically. The device memory can store up to 192 results for each of the two users. Select memory

for users  and  is made before starting the measurement when the device is switched off by pressing the button  as many times as necessary to select the correct memory.


The following phases of the measurement are shown in the illustrations below.

Press button "ON/OFF".

All symbols in the display switch on for 1 second.

The screen in the DIA column should display "0" (zero), which means that the blood pressure monitor is ready to start measuring.

HOFFEN

EN

The measurement will start automatically by inflating the cuff. The cuff is pumped up to 180-190 mmHg as standard. If a higher pressure is required to take the measurement, the device increases the pressure on pressing

the button . Pressing this button allows to choose between two pumping pressure values: 260 and 280 mmHg

The maximum pressure to which the cuff can be inflated is 280 mmHg. We recommend that the cuff inflation value should be a maximum of 20-30 mmHg higher than expected measurement value.

A rising pressure value will appear on the display until the maximum cuff inflation value is reached.

CAUTION! *If you are experiencing major discomfort during cuff inflation, or if the cuff inflation does not end when the set top value is reached, immediately release and remove the cuff and switch off the device with the "On/Off" button.*

When you start deflation of the cuff, a symbol  will appear indicating that a pulse has been detected.

When the measurement is completed correctly, the screen will display the values - from top: systolic (SYS - in this case 128 mmHg), diastolic (DIA - in this case 78 mmHg) and pulse (PULSE - in this case 72 times per minute) The result of the correct measurement will be stored in the memory determined for the user.

If the user does not sit still during the measurement (moves the entire body or the hand on which the cuff is placed), two symbols may appear on the screen:  or 

. This means that the measurement is incorrect and must be repeated. We recommend resting for at least 3 minutes and re-measuring.

If the measurement shows an irregular heartbeat, the symbol  will be displayed. Irregular heartbeat is a condition where the heartbeat is 25% slower or faster than the average heartbeat measured during systolic and diastolic blood pressure measurement. We recommend resting for at least 3 minutes and re-measuring. If the symbol appears again, it is necessary to contact your doctor.

Remove the cuff after the measurement is taken. Switch off the unit by pressing the "ON/OFF" button. The device can also be left standing and will turn off automatically after 2 minutes.

VIEWING THE RESULTS STORED IN THE DEVICE MEMORY

The memory of the device can store 192 measurement results for each of the 2 users.

You can look through the stored measurements. To do this, press the button  while the device is switched off and use it to select whether to display the results for the user

 or  for the user.

Pressing the button  while the device is switched off will enable viewing the saved results (of the current user).


Then use the  and  buttons to view the stored data.

Switch off the unit by pressing the "ON/OFF" button when viewing is complete. When the memory is full (192 measurements taken for one user), each subsequent result will be stored starting from position 0 (zero) and the oldest result shall be deleted at the same time. The list will move by one position.

ADDITIONAL SETTINGS

The user has the possibility of additional setting of some parameters of the device. It can select the measurement units in which the result is to be displayed - mmHg and kPa (1kPa = 7.5 mmHg), activate or deactivate voice messages and clear the measurement memory. To perform changes, you need to run an additional device menu. To do this on a switched off device press and hold the button for 3 seconds . Pressing it again will toggle the menu, as shown in the following figure.

HOFFEN


The symbol  around an item in the above and subsequent illustrations indicates that it is blinking.

WIPING OUT THE DATA

It is not possible to delete one selected measurement.

In order to clear all measurements for the current user (previously set), select  and press the button  in the additional settings menu. The confirmation of erasing the measurement memory will be the change of the flashing value corresponding to the number of stored results (in this case the number 3) to the flashing number 0 (zero), as shown in the following figure.


CHANGE OF UNITS OF MEASUREMENT

To change measurement units from mmHg to kPa or vice versa, select  mmHg or kPa (depending on which unit was previously selected) and change the selection with the button  in the additional settings menu. The selection is confirmed by pressing

HOFFEN


EN

the button  again or turning off the unit by pressing the "ON/OFF" button. The change is shown in the illustration below.


TURNING THE SOUND ON AND OFF

In order to disable or activate the sound (messages) in the device, enter the additional menu select with  the icon  or button , depending on the previous settings. When the icon starts flashing, change the selection with . The selection is confirmed by pressing the button  again or turning off the unit by pressing the "ON/OFF" button. The change is shown in the illustration below.

POWER SUPPLY FROM AN EXTERNAL POWER SOURCE


 In order to be able to use an external power supply, you must remove the entire set of batteries from the battery compartment.

The blood pressure monitor can be powered not only by battery power, but also by an external power supply via the microUSB socket. Simply connect the charger to your phone or smartphone, powerbank or computer with 5V – 500mA current at its output (not included), through cable with a plug compatible with these energy sources and a microUSB plug-in.

WHO (WORLD HEALTH ORGANIZATION) BLOOD PRESSURE STANDARDS

The chart below presents the classification of blood pressure according to WHO.

9. CLEANING AND MAINTENANCE

Cleaning should only be performed when the unit is switched off.

Switch off the unit by pressing the "ON/OFF" button before cleaning. The outside of the unit and the rubber air tube may only be cleaned with a dry cloth or paper towel. Never immerse the unit in water or wash it under running water. In case of heavy soiling, you can wipe the outer parts of the unit with a damp, well wrung cloth and then immediately wipe it dry with a paper towel or a dry cloth. It is not allowed to use corrosive liquids, petrol, solvents, etc. for cleaning the device.

The cuff can be cleaned with a damp, heavily wrung cloth and a small amount of ordinary soap (if the cuff is heavily soiled). After cleaning, leave the cuff to dry.

10. TROUBLESHOOTING

In case of any problems with the device, please refer to the following tips.

<p>The device does not switch on.</p>	<p>Check that the battery is inserted correctly - that the polarity is maintained. Check the batteries. If they does not work properly, replace the whole set with a new one.</p>
<p>The following symbol  appears on the display.</p>	<p>Batteries low, replace the whole set with a new one,</p>
<p>The device is working but cannot inflate the cuff and the display shows the symbol </p>	<p>Check that the air hose is securely plugged into the outlet on the device, or that there is no dirt between the air tube plug and the outlet on the unit that could cause the air to leak. The cuff is applied to the clothing and not directly to the skin. The cuff is incorrectly attached or not sufficiently tightened with velcro. The user moved the arm during the measurement. The inflation pressure is not higher than the systolic pressure. Check the air tube or cuff for damage.</p>
<p>The following symbol  appears on the display.</p>	<p>Cuff is not connected or is not properly connected to the device. Air leakage at the socket of the device, airway tube or cuff.</p>
<p>The following symbol  appears on the display.</p>	<p>Measurement error</p>
<p>The following symbol  appears on the display.</p>	<p>Pressure higher than 280 mmHg - consult a doctor!</p>


<p>The following symbol</p>  <p>appears on the display.</p>	<p>Appears when you want the preview of the measurement memory for the user in question, but there is no record yet.</p>
<p>The reading result is too low (or too high).</p>	<p>The cuff is applied to the clothing and not directly to the skin. The cuff is incorrectly attached or not sufficiently tightened with Velcro. The user moved the arm during the measurement.</p>
<p>Cuff pressure's dropping too fast.</p>	<p>Cuff's too tight.</p>
<p>Other errors</p>	<p>Switch off the unit by pressing the "ON/OFF" button and switch it on again. Replace the entire set of batteries with new, efficient ones.</p>
<p>If the problem persists, contact the manufacturer's service department.</p>	

11. STORAGE AND TRANSPORT

STORAGE OF EQUIPMENT

If you do not use the device for a short period of time, store it in the supplied case after turning it off and cleaning it if necessary. Try not to twist and fold the air tube and cuff too tightly.

If you are planning not to use the device for a long time clean it, disconnect the air tube and put it in a slightly rolled-up cuff (as shown in the illustration), remove the battery pack, store it in the supplied case in a dry room away from heat and moisture sources. It is best to place the device in its original packaging, away from children and pets.

TRANSPORT OF THE DEVICE

When handling, only carry the unit after

it has been switched off, holding it securely by the housing and the cuff.

It is recommended to transport the device only when the following conditions have been met: Clean the device, remove the battery pack and disconnect the air tube from the device. Put the device into the case and preferably into the original packaging. Do not expose the unit to excessive vibration, shock, heat, humidity or extreme temperatures during transportation.

12. WARRANTY AND SERVICE

The product is covered by a 3 year warranty from the date of sale of the product. Batteries are not covered by the warranty - it is a consumable material. If you need warranty repair, please contact the manufacturer's hotline. The product returned for repair should be complete and preferably in its original packaging.

13. DECLARATION OF CONFORMITY WITH EUROPEAN STANDARDS

 This device meets the requirements of the European Commission Directive 93/42/EEC on Medical Devices. This blood pressure monitor is designed in accordance with European Standard EN 1060, Non-invasive Sphygmomanometers, Part 1, General Requirements, and Part 3.

14. MANUFACTURER AND IMPORTER:

MANUFACTURER:

Cofoe Medical Technology Co., Ltd.
816 Zhenhua Road, Yuhua District
410000 Changsha, Hunan, PEOPLE'S REPUBLIC OF CHINA
SUNGO Europe B.V.
Olympisch Stadion 24, 1076DE Amsterdam

IMPORTER:

mPTech Sp. z o.o.
Nowogrodzka 31 street
00-511 Warsaw
Poland
Helpline: (+48 71) 71 77 400

Manufactured in China

LOT 612736

Date of issue of the manual: 2022.06

15. ELECTROMAGNETIC COMPATIBILITY INSTRUCTIONS

- Automatic brachial blood pressure monitor meets the requirement of electromagnetic compatibility in IEC60601-1-2.
- The Automatic brachial blood pressure monitor (KF-65A) is suitable for home healthcare environments.
- The user needs to install and use according to electromagnetism compatibility information which is attached with it.
- Guidance and manufacture's declaration stated in the appendix A.

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this Automatic brachial blood pressure monitor (KF-65A) adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, the equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment may influence performance, it should be used no closer than 30 cm (12 inches) to any part of Automatic brachial blood pressure monitor (KF-65A). Otherwise, degradation of the performance of this equipment could result.

Appendix A. EMC Information

Guidance and manufacturer's declaration – electromagnetic emissions

Emission test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not application
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not application

Guidance and manufacturer's declaration – electromagnetic immunity

Immunity test	IEC 60601-1-2 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air
Electrical fast transient/burst IEC 61000-4-4	Not application	Not application
Surge IEC 61000-4-5	Not application	Not application
Voltage dips, short interruptions and voltage	Not application	Not application

HOFFEN

EN

variations on power supply input lines IEC 61000-4-11		
Power frequency magnetic field IEC 61000-4-8	30 A/m 50/60 Hz	30 A/m 50/60 Hz
Conducted RF IEC 61000-4-6	Not application	Not application
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz
NOTE: U_T is the AC mains voltage prior to application test level.		

Guidance and manufacturer's declaration – electromagnetic immunity

Radiated RF IEC 61000-4-3

(Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)

Test frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1,8	0,3	27
450	430-470	GMRS 460, FRS 460	FM +/- 5 kHz Deviation 1 kHz sine	2	0,3	28
710	704-787	LTE band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE pasmo 5	Pulse modulation 18 Hz	2	0,3	28
870						
930						
1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE band 1, 3, 4, 25, UMTS	Pulse modulation 217 Hz	2	0,3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE pasmo 7	Pulse modulation 217 Hz	2	0,3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9
5500						
5785						